

Vaporized Hydrogen Peroxide (VHP®) Gaseous Decontamination : 'GREEN' Technology for the Highest Level of Microbial Control within a Pharmaceutical Facility

Peter Harris, B & V TESTING
Larry Zanko, STERIS Corporation

ISPE Boston Area Chapter
February 23, 2010



ENGINEERING PHARMACEUTICAL INNOVATION



Agenda – Part 1

- Overview of the VHP Process: efficacy, safety, regulatory landscape, environmental properties and material compatibility
- VHP Technology and Equipment
- VHP Biodecontamination Field Applications
- VHP Services and Project Planning
- Case Studies

ENGINEERING PHARMACEUTICAL INNOVATION



Why Gaseous Decontamination?

- Contamination remediation
- Pre-occupancy new or renovated facility
- Preventative periodic bio-burden reduction
- Product/population change
- Equipment transfer
- Facility decommissioning
- Equipment maintenance, i.e. BSCs, HEPA housings, etc.

ENGINEERING PHARMACEUTICAL INNOVATION




VHP: Background and History

- VHP process developed mid-1980s, utilizing patented closed-loop, low concentration “dry” process
- VHP process patent issued in early 1990s
- Over 1,200 VHP systems in use world-wide in pharma, med device, LAR, bio-containment fields
- 15+ years of validated use in pharmaceutical manufacturing
- Large scale VHP facility applications applied for remediation during 2001 anthrax attacks
- 2007 VHP contract service offering available to industry

ENGINEERING PHARMACEUTICAL INNOVATION




VHP: the Process



ENGINEERING PHARMACEUTICAL INNOVATION



“Typical” VHP Biodecontamination Cycle



ENGINEERING PHARMACEUTICAL INNOVATION



Variables Affecting Efficacy

- Temperature and humidity affect how much hydrogen peroxide (HP) can be generated in the gas state, start Rh typically 30-50%
- Concentration—typically 150-400 ppm room applications
 - (injection rate/ air flow) * wt. % of HP
- Saturation
 - Inject as much as possible but below dew-point
 - Humidity is good for microbial kill
- Distribution—may be facilitated with fans
- Materials which can reduce VHP concentration, cellulosic material (i.e. cardboard, paper), galvanized steel, standing water

ENGINEERING PHARMACEUTICAL INNOVATION



Why H₂O₂ Vapor? Sporicidal at Low Concentrations

- Bacterial Spores - Most Resistant Organism to VHP
- Highly sporicidal even as low as 0.1 mg/l
- Broad kill spectrum
- 35% H₂O₂ Registered EPA sterilant
- As per EPA label short “contact” times, **six log reduction: 30-90 minutes @ 250-400 ppm**




Geobacillus stearothermophilus
Cross section

ENGINEERING PHARMACEUTICAL INNOVATION



Resistance to VHP for Biological Organism Classes



ENGINEERING PHARMACEUTICAL INNOVATION



Why H₂O₂ Vapor? Efficacy Testing and Monitoring

- Biological Indicators - Most VHP Resistant *geobacillus stearothermophilus* typically 4-6 log
- Chemical Indicators – Qualitative instant feedback
- Real-time monitoring with electro-chemical sensors




ENGINEERING PHARMACEUTICAL INNOVATION



VHP Kill Matrix

Geobacillus stearothermophilus spores inoculated on Stainless Steel Coupons at 30°C



ENGINEERING PHARMACEUTICAL INNOVATION



Bacterial Spores Evaluated for Their Resistance to Hydrogen Peroxide Gas

- *Bacillus anthracis*
- *Bacillus cereus*
- *Bacillus circulans*
- *Bacillus pumilus*
- *Geobacillus stearothermophilus*
- *Bacillus subtilis*

- *Clostridium botulinum*
- *Clostridium sporogenes*
- *Clostridium tetani*
- *Clostridium difficile*

ENGINEERING PHARMACEUTICAL INNOVATION



Mycobacteria Evaluated for Their Resistance to Hydrogen Peroxide Gas

- *Mycobacterium smegmatis*
- *Mycobacterium terrae*
- *Mycobacterium bovis*
- *Mycobacterium tuberculosis*
- *Nocardia lactamdurans*

ENGINEERING PHARMACEUTICAL INNOVATION



Non-enveloped, Non-lipid Viruses Evaluated for Their Resistance to Hydrogen Peroxide Gas

- *Parvoviridae* (*Feline and Canine parvovirus*)
- *Picornaviridae* (*Polio Type 1, Swine Vesicular, Rhinovirus 14*)
- *Reoviridae* (*Blue Tongue, Avian reovirus*)
- *Caliciviridae* (*Vesicular exanthema*)

ENGINEERING PHARMACEUTICAL INNOVATION



Gram-negative Vegetative Bacteria Evaluated for Their Resistance to Hydrogen Peroxide Gas

- *Burkholdia cepacia*
- *Pseudomonas aeruginosa*
- *Serratia marcesens*
- *Escherichia coli*
- *Proteus vulgaris*
- *Salmonella choleraesuis*

ENGINEERING PHARMACEUTICAL INNOVATION



Fungi, Molds, and Yeasts Evaluated for Their Resistance to Hydrogen Peroxide Gas

- *Aspergillus niger*
 - *Aspergillus terrus*
 - *Candida parapsilosis*
 - *Rhodotorula glutinis*
 - *Fusarium oxysporum*
 - *Penicillium chrysogenum*
- *Candida parapsilosis*
 - *Saccharomyces cerevisiae*
 - *Rhodotorula glutinis*

ENGINEERING PHARMACEUTICAL INNOVATION



Large Non-enveloped Viruses Evaluated for Their Resistance to Hydrogen Peroxide Gas

- *Adenovirus (Adenovirus 2)*
- *Poxviridae (Vaccinia)*

ENGINEERING PHARMACEUTICAL INNOVATION



Gram-positive Bacteria Evaluated for Their Resistance to Hydrogen Peroxide Gas

- *Enterococcus faecium*
- *Enterococcus faecalis*
- *Staphylococcus aureus* (MRSA)
- *Lactobacillus casei*
- *Listeria monocytogenes*
- *Legionella pneumophila*

ENGINEERING PHARMACEUTICAL INNOVATION




Enveloped, Lipid Viruses Evaluated for Their Resistance to Hydrogen Peroxide Gas

- *Orthomyxoviridae* (*Avian Influenza*)
- *Paramyxoviridae* (*New Castle*)
- *Herpesviridae* (*Pseudorabies, Herpes Simplex*)
- *Rhaboviridae* (*Vesicular stomatitis*)
- *Toga/Flaviviridae* (*Hog cholerae, BVD*)

ENGINEERING PHARMACEUTICAL INNOVATION



Why H₂O₂ Vapor? Safety Profile



Hydrogen
Peroxide

Skin/eye irritant
PEL 1.0 ppm
IDLH 75 ppm

Chlorine Dioxide

Severe irritant
PEL 0.1 ppm
STEL 0.3 ppm
IDLH 5.0 ppm

Formaldehyde


Human carcinogen
PEL 0.75 ppm
STEL 2 ppm
IDLH 10 ppm

PEL - Permissible Exposure Limit (8 hours)
STEL - Short-term Exposure Limit (15 min.)
IDLH - Immediately Dangerous to Life or Health

ENGINEERING PHARMACEUTICAL INNOVATION



Why H₂O₂ Vapor? Safety Profile




Typical Decontamination Concentrations:


VHP	Chlorine Dioxide	Formaldehyde
150-400 ppm	500-1,500 ppm	8,000-10,000 ppm

Why H₂O₂ Vapor? Safety Profile, VHP Degradation: real-time monitoring

VHP Concentration



Why H₂O₂ Vapor? Safety Profile



When considering VHP typical use concentrations, degradation and migratory properties and OSHA exposure limits, VHP offers less exposure risk to personnel and animals where small leaks occur.

Why H₂O₂ Vapor? A Green Solution

Non-toxic byproducts

No post-decon wipe down



EPA: “no risks to the environment are expected from use of pesticide products containing hydrogen peroxide because 1) the substance readily decomposes to water and oxygen gas, leaving no residue; 2) it is effective at low concentrations where no toxic effects are expected” www.epa.gov

Material Compatibility*

Metals

Aluminum Excellent
Anodized Aluminum Good
Brass Good
Copper Good.
Stainless Steel – all grades Excellent
Steel Good
Titanium Excellent

Stainless Steel – all grades Excellent
Steel Good
Titanium Excellent

Plastics

ABS Excellent
Aflas Excellent
CPVC (chlorinated polyvinylchloride) Excellent
Kel –F Excellent
Nylon Fair
PMMA Excellent
Polyethersulfone (PES) Excellent
PolyethyleneTerephthalate (PET) Good

Polyacetal (Delrin) Excellent
Polycarbonate Excellent
Polyetherimide (Ultem) Excellent
Polymethylepentene Excellent
Polyphenylene oxide Excellent
Polyetherketone (PEEK) Excellent
Polyethylene (HDPE, LDPE UHMWPE) Excellent

Elastomers

Buna N Fair
Butyl Rubber Good
Chem Raz Good
EPDM Fair
Viton Excellent

Hypalon Excellent
Polyurethane Excellent
Silicone Rubbers Excellent
Kelrez Excellent

*Compatibility is defined as the materials ability to undergo exposure to VHP with no significant changes in physical, or chemical properties (i.e. no changes in strength, flexibility, chemical composition, color etc.).

ENGINEERING PHARMACEUTICAL INNOVATION



Why H₂O₂ Vapor?

Material and Component Compatibility Pharma Project Pilot Field Test


- laptop computer (turned on)
- LCD monitor (turned off)
- telephone (on)
- electronic scales (on)
- various electronic sensors
- stainless steel clamps, fittings, and connectors
- rubber grommets and washers
- rubber hoses
- plastic funnel and containers



ENGINEERING PHARMACEUTICAL INNOVATION



Why H₂O₂ Vapor? EPA-registered, FIFRA compliant



In compliance with Federal
Insecticide, Fungicide and
Rodenticide Act, all anti-microbial
agents must be EPA registered.
35% hydrogen peroxide
EPA reg. no. 58779-4



Why VHP? A Summary of Advantages

- Rapid decontamination at ambient temperatures and low concentrations
- Strong history of use and efficacy data, easily validated with BIs
- Strong comparative safety profile—detectable well below IDLH, no “trapped gases”
- Non-toxic by-products—a “green” solution
- No lengthy aeration
- No residue
- Strong material and component compatibility profile
- EPA Registered—FIFRA compliant

VHP Pharmaceutical Facility Applications

Rapidly Increasing in Popularity, Common Applications Include:

- Aseptic Manufacturing Rooms and Pilot Production Rooms: pre-occupancy new or renovated facility, remediation of know contaminant, periodic preventative
- Tissue Culture Rooms, Cold Rooms, Warm Rooms
- Lab Animal Research Procedure, Cage Change and Equipment Transfer Rooms
- Biological Safety Laboratories Level 3 (BSL-3)
- Biological Safety Cabinets, Incubators, Enclosures

ENGINEERING PHARMACEUTICAL INNOVATION



VHP Generation Technology



- Mobile VHP
- Installed VHP
- Facility integrated VHP
- EPA-registered 35% HP
- For Room Applications with ARD use, for 3,500 ft³ room, typically 1 machine will achieve six-log kill in 2 hours or 7,000 ft³ in 4 hours

ENGINEERING PHARMACEUTICAL INNOVATION



Methods of VHP Injection


- External via integrated ports
- External via installed Plexiglas panels
- In room, solo or daisy-chained
- Via AHUs (can include ductwork decontamination)



ENGINEERING PHARMACEUTICAL INNOVATION




VHP Room Decontamination Typical Room Application



ENGINEERING PHARMACEUTICAL INNOVATION



Flow Diagram



ENGINEERING PHARMACEUTICAL INNOVATION



VHP Service Project Planning

- Fumigation management plan (FMP) project document
- Define purpose and scope of decontamination
- Identify the team players including stakeholders and support personnel (area managers and PIs, EH &S, facilities engineering, security, etc.) and assign responsibilities
- Review site schematics, perform site visit, evaluate HVAC and electrical capabilities
- Review personnel/authority notification, site control and security, site signage

ENGINEERING PHARMACEUTICAL INNOVATION



VHP Project Planning contd.

- Review area preparation: pre-cleaning, material/equipment transfer, HVAC control/support, smoke detector disengagement responsibilities, sealing of space
- Establish safety buffer zone, project safety plan and external monitoring plan
- Review post decontamination area clearance and equipment/material retrieval procedures
- Define acceptance criteria, establish BI quantity/mapping where applicable
- Establish the final project schedule, task sequencing and responsibilities

ENGINEERING PHARMACEUTICAL INNOVATION



VHP Project Execution

- Upon arrival to job site equipment, material and personnel transfer commences in accordance with site requirements (gowning, equipment transfer protocols)
- Target area is prepared for decontamination: BI, CI, VHP generator, fan, emergency signage placement
- Smoke detectors disengaged/ HVAC isolated
- Final pre-go walk-thru/assessment with client to confirm area readiness
- Upon client authorization/clearance VHP injection commences
- Real-time internal and external monitoring of VHP
- Following injection commensurate with project requirements, aeration commences to bring VHP levels in area to below PEL 1.0 ppm
- Following area clearance, equipment, and BI/CI retrieval
- Upon completion of independent third-party BI analysis, report issued

ENGINEERING PHARMACEUTICAL INNOVATION



Case Study 1: Pharmaceutical Manufacturing Facility

- Project scope: 240,000 ft³ manufacturing space including bioreactor rooms
- Issue: bacterial remediation emergency
- Remote (AHU) and terminal ceiling HEPA filtration
- AHU and local VHP injection




VHP Case Study 1 continued

- 65 rooms, 48 hours
- 56 liters of Vaprox 35%
- 165 biological indicator locations
- 95 Chemical Indicators
- VHP injection 12 g/min
- VHP average concentration 242 ppm



VHP Case Study 1 continued



ENGINEERING PHARMACEUTICAL INNOVATION



VHP Case Study 1: Results

- All CIs changed color
- 96% BIs ≥ 6 log kill, with remaining ≥ 5 log kill
- External monitoring VHP concentration < 1.0 ppm throughout safety perimeter



ENGINEERING PHARMACEUTICAL INNOVATION



Case Study 2: Pharmaceutical Manufacturing Facility

- Project Scope: 84,000 ft³ pilot production facility
- Issue: Preventative during shutdown
- Remote (AHU) and terminal ceiling HEPA filtration
- Local VHP injection



ENGINEERING PHARMACEUTICAL INNOVATION



VHP Case Study 2 continued


- 28 rooms, 24 hours
- 22 liters of Vaprox 35%
- 65 biological indicator locations
- 28 Chemical Indicators
- VHP injection 12 g/min
- VHP average concentration 282 ppm



ENGINEERING PHARMACEUTICAL INNOVATION



VHP Case Study 2 continued



ENGINEERING PHARMACEUTICAL INNOVATION



VHP Case Study 2 continued

- All CIs changed color
- Six-log kill all 65 BIs
- External monitoring VHP concentration < 1.0 ppm throughout safety perimeter



ENGINEERING PHARMACEUTICAL INNOVATION



Questions and Project Discussion

Contacts:

Peter Harris peter@bandvtesting.com
Nick Flynn nick.flynn@bandvtesting.com,
800-851-9081
www.bandvtesting.com
www.steris.com

ENGINEERING PHARMACEUTICAL INNOVATION



Agenda – Part 2

- H_2O_2 - Vapor vs. Condensate
- Portable or Modular?
- Installation of Modular VHP Systems
- Modular VHP Case Studies
- Final Thoughts

ENGINEERING PHARMACEUTICAL INNOVATION



Vapor or Condensate?

Boiling Points:
 H_2O 100°C
 H_2O_2 150°C




...if you can see it,
it's not a vapor

ENGINEERING PHARMACEUTICAL INNOVATION



Vapor (300ppm) vs. Condensate (700,000ppm)



Effect of
improper
application

ENGINEERING PHARMACEUTICAL INNOVATION



VHP Flexibility



- Maximum flexibility – multiple enclosure types
- Rooms, product or package handling equipment
- Integrations with Chambers, Washers and Autoclaves



ENGINEERING PHARMACEUTICAL INNOVATION



Which VHP System?

Portable



- ✓ Spaces not yet defined
- ✓ Uses in different buildings
- ✓ Typically less than 10,000ft³
- ✓ Cycle time not a constraint
- ✓ Use of fans not an issue
- ✓ Less frequent use

Modular



- ✓ Large and small spaces up to ~80,000 ft³
- ✓ Same enclosures repeatedly
- ✓ Frequent use (chamber)
- ✓ Short cycle times
- ✓ Automated sequenced decontamination of multiple rooms

ENGINEERING PHARMACEUTICAL INNOVATION



Why Modular?

Keep Equipment Outside Space

- Save space within room / Pass Through
- Avoid cross contamination
- Keep maintenance activities outside

No Set Up

- Decon at the Press of a Button
- Run Sequential Decons via BMS*
- Reduced Handling of Peroxide Excellent Distribution

Cost

- Less expensive than multiple portables
- Save on labor
- The easier to use – the more frequent the use -the cleaner the space


ENGINEERING PHARMACEUTICAL INNOVATION



*BMS = Building Management System

VHP Modular Integration Schematic

Single Pass



VHP Exhaust from Room

ENGINEERING PHARMACEUTICAL INNOVATION



VHP Emissions

- No EPA Limit
- Uniform Fire Code
 - ½ IDLH emergency conditions for “toxic” gases
 - Local regulations ? Check.

Exhaust is rarely if ever an issue.

Why?

- Dilution
- Break down (galvanized ductwork)
- Easy to install catalyst



ENGINEERING PHARMACEUTICAL INNOVATION



There are **2** options for integration...

Airflow	Single Pass	Recirculating
Typical Applications	BSL, Lab Animal	Clean Room, RABS
Dehumidification	HVAC	HVAC
VHP Supply	Dedicated each room	HVAC
Pressure control	Dedicated blower	HVAC
Aeration	HVAC	HVAC




T4

ENGINEERING PHARMACEUTICAL INNOVATION



Single Pass System




ENGINEERING PHARMACEUTICAL INNOVATION



Single Pass System


- VHP Pipes
 - Made of CPVC or PP
 - Insulated (supply only) but not traced
- Cycle Phases
 - Dehumidification via HVAC
 - Decontamination phase
 - HVAC is stopped
 - Leaktight dampers are closed
 - Gas flows into and from the rooms at specified rates via pre-set butterfly valves
 - Aeration Phase
 - HVAC restarted (VHP - exhausted to outside)
- Pressure Control (Positive or Negative)
 - Possible during decontamination cycles
 - On small and/or leaktight rooms
 - By variable speed exhaust fan



ENGINEERING PHARMACEUTICAL INNOVATION



Recirculating System (Clean Room)



Recirculating System

- VHP Pipework
 - Made of CPVC or PP.
 - Section to central supply insulated but not traced
- Ductwork / HEPAs
 - Must be air-tight
 - Materials - PVC coated, galvanized, Aluminum , Stainless Steel but no bare copper
 - VHP will penetrate and decontaminate HEPAs
- Cycle Phases:
 - Dehumidification/ Injection Phase
 - Runs closed loop only , with no fresh air.
 - Recirculation allows an even distribution of the gas concentration in all the rooms .
 - Cooling and heating of the HVAC will be stopped .
 - Aeration: the maximum of fresh air is admitted during the Aeration Phase

Large Room Biopharmaceutical Fermentation Suite

Volume: 32,000ft³ (900m³)

Ceiling height: 28ft (8,5m)

Single pass, No fans

6-log reduction

Cycle time 6 hours

New construction

Cycle Phase	Time min.	Airflow	Injection g/min
Dehumidification	30	6 A.E./ hour	-
Condition	30	120 cfm	96
Decontamination	90	120 cfm	60
Aeration	210	40 A.E./ hour	-




ENGINEERING PHARMACEUTICAL INNOVATION



A.E. = Air Exchange

RABS – (Restricted Access Barrier Systems)

Modular VHP systems can rapidly decontaminate both Active & Passive RABS and the surrounding rooms.



Source: Pharmaceutical International

ENGINEERING PHARMACEUTICAL INNOVATION



Pass-Through Chambers / Transfer Hatches / Material Air Locks



Shown
above

Enclosure v Volume ft ³	Enclosure surface material	Injection rate Condition g/min.	Injection rate Decon g/min.	Decon time Min 6 log	Decon airflow ft ³ /min	PPM	Aeration airflow ft ³ /min	Total Cycle Time min.
460 (6x8x9.5'L)	Stainless	32	23	12	120	1000	765	45
175 (4x6x7'L)	Epoxy paint	12	9	8	40	950	1750	30

ENGINEERING PHARMACEUTICAL INNOVATION



HEPA Filter Decontamination



ENGINEERING PHARMACEUTICAL INNOVATION




BSL Lab & BSC Decon Single Pass


**Simultaneous Decon
of Primary Containment**

A2 type biosafety cabinets can be decontaminated together with the room

- Exhaust dampers above cabinets are closed
- Cabinet blowers left on



ENGINEERING PHARMACEUTICAL INNOVATION



Automated Sequential Zone BSL-3 Lab Decontamination Single Pass

The floor plan illustrates the layout of a BSL-3 laboratory, divided into seven distinct zones for sequential decontamination. Each zone is color-coded and labeled:

- Zone 1:** Light blue, located at the top left.
- Zone 2:** Yellow, located below Zone 1.
- Zone 3:** Purple, located at the top right.
- Zone 4:** Pink, located below Zone 3.
- Zone 5:** Light green, located at the bottom left.
- Zone 6:** Green, located below Zone 5.
- Zone 7:** Dark pink, located at the bottom right.

Key features and labels on the floor plan include:

- Emergency Button:** Located at the top left, with text: "EMERGENCY BUTTON FOR INITIATION OF DECONTAMINATION SYSTEM. SIGNALS TO FILL MISTERS, AIRWAYS, FLOWERS, AND MISTERS. WILL ALSO PROVIDE OVERSIGHT FOR FILL MISTERS TO BE USED. THIS FILLING SYSTEM IS USED TO FILL MISTERS TO REAL. EMERGENCY DECONTAMINATION. USE OF THIS SYSTEM MUST BE INITIATED AT 10' ABOVE FLOOR TO CENTER OF BUTTON. TYPICAL OF 8 LOCATIONS."
- Alarm Horn and Smoke Detector:** Located at the bottom left, with text: "ALARM HORN AND SMOKE DETECTOR. EMERGENCY DECONTAMINATION SYSTEM. SIGNALS TO FILL MISTERS, AIRWAYS, FLOWERS, AND MISTERS. WILL ALSO PROVIDE OVERSIGHT FOR FILL MISTERS TO BE USED. THIS FILLING SYSTEM IS USED TO FILL MISTERS TO REAL. EMERGENCY DECONTAMINATION. USE OF THIS SYSTEM MUST BE INITIATED AT 10' ABOVE FLOOR TO CENTER OF BUTTON. TYPICAL OF 8 LOCATIONS."
- Manhole:** Located at the bottom right, with text: "MANHOLE. SIGNALS AND FILLING TO FILL MISTERS. TYPICAL OF 8 LOCATIONS."
- Zone 7:** Located at the bottom right, with text: "ZONE 7. FILL MISTERS TO BE INITIATED."

The photograph on the right shows the physical hardware of the automated decontamination system, featuring several vertical pipes with valves and sensors, labeled "MANHOLE" and "FILL MISTERS".

Zone

- 1
- 2
- 3
- 4
- 5
- 6
- 7

ENGINEERING PHARMACEUTICAL INNOVATION

ISPE

Task List

Vendor	General Contractor
<ul style="list-style-type: none"> ✓Parts <ul style="list-style-type: none"> VHP Generator Dehumidifier (Munters) Bulk Fill apparatus Extra Controllers (remote) Sensors/ Monitors ✓Consulting / Labor <ul style="list-style-type: none"> System Selection / Layout VHP Air Supply Balance Cycle Development Interface with BMS IQ/OQ and commissioning VHP Unit Installation 	<ul style="list-style-type: none"> ✓Parts <ul style="list-style-type: none"> Ductwork (piping) Dampers (airtight) Booster Fan (if needed) ✓Labor <ul style="list-style-type: none"> BMS Integration Ductwork Installation

ENGINEERING PHARMACEUTICAL INNOVATION



Selected Installations

Pharma / Animal Health

- GSK
- Intervet
- Fresenius
- Sanofi
- Merial
- Alcon
- Pfizer
- Fort Dodge

Public Health Labs

- Indiana
- New Jersey
- West Virginia

Others

- Tripler- US Army
- Univ. Nebraska
- Lawrence Livermore
- National Labs
- INRS



PHARMACEUTICAL INNOVATION



Why VHP?

➤ Environmentally Safe

- ✓ Excellent Material Compatibility, Even with Sensitive Electronics, and Low Toxicity

➤ Residue-Free

- ✓ Quickly Breaks Down into Water Vapor and Oxygen.

➤ Ideal for Cleanrooms

- ✓ Integrated as a user-friendly Utility

➤ Use Registered with EPA



STERIS®

ENGINEERING PHARMACEUTICAL INNOVATION



Why VHP?

✓ Consistency & Distribution

- Wet surfaces / minimal contact times -not an issue
- Passes through HEPA filters
- Decontaminates biosafety cabinets and HEPAs during room decon
- Rapidly kills airborne and surface microbes

✓ Labor

- Minimal labor required
- Hundreds of validation applications



Reduce



ENGINEERING PHARMACEUTICAL INNOVATION



Acknowledgements

Claire Fritz and John Klostermyer
STERIS Corporation

ENGINEERING PHARMACEUTICAL INNOVATION



Thank you!

Contacts:

Peter Harris peter@bandvtesting.com

Nick Flynn nick.flynn@bandvtesting.com

800-851-9081

www.bandvtesting.com

Contacts:

Larry Zanko larry_zanko@steris.com

Leena Asplund leena_asplund@steris.com

800-548-4873

www.steris.com

ENGINEERING PHARMACEUTICAL INNOVATION

