



# INTEGRATED VPHP DECONTAMINATION SYSTEMS *THE EMERGING UTILITY*

John Klostermyer, Ph.D.  
ISPE Product Show  
Track 3 Session 3  
September 26, 2018

## AGENDA

1. BACKGROUND INFORMATION
2. BENEFITS
3. PROCESS
4. BEST MANAGEMENT PRACTICES
5. VALIDATION / SAFETY
6. STAKEHOLDERS
7. FAQs

### TECHNICAL—FACILITIES AND EQUIPMENT

## INTEGRATED VPHP DECONTAMINATION SYSTEMS: THE EMERGING UTILITY

John Klostermyer, Ph.D., Albert Garcia, and Don E. Morgan

Integrated VPHP systems offer a versatile, automated, specialized process for decontaminating surfaces, isolators, RABS, chambers, and pass-through hoods.

The first portable vapor-phase hydrogen peroxide (VPHP) generators developed in the early 1990s were designed to dry, decontaminate, and sanitize target enclosures efficiently while controlling pressure. Today, these generators, which typically use a closed-loop airflow pattern, are used predominantly on isolators, small rooms, and material air locks. As users sought to enhance process automation, increase decontamination rates, shorten cycle times, and standardize data capture, integrated VPHP generators were developed. These could be tied in with the heating, ventilation, and air conditioning (HVAC) system and controlled and monitored by the building management system (BMS).

Integrated generators are designed to rapidly absorb hydrogen peroxide precursors and consistently deliver regulated quantities of VPHP to the target enclosure. Unlike portable VPHP generators, integrated units work with an array of internal air-handling components that deliver, distribute, and purge VPHP from the target enclosure. Collectively, this is called an integrated VPHP system. In some facilities, use of an integrated VPHP system is common for other duties, such as controlled air.

Integrated VPHP installations provide the well-known benefits of hydrogen peroxide vapor decontamination while enhancing process repeatability and decreasing labor and total costs. Other key benefits include:

- Rapid cycle surface decontamination
- No residues
- Excellent material compatibility
- Lower toxicity than other gaseous treatments
- Highly automated, high-yield, consistent, and continuous operation
- Very low operating cost, no manual labor effort, or time
- Systems installed and maintained outside the clean area
- No cross-contamination via equipment or personnel moving between target enclosures
- Single unit can be configured via a manifold to decontaminate multiple enclosures

**BACKGROUND**  
Initial VPHP applications were focused predominantly on pharmaceutical sterilization and aseptic production applications by purging isolators with portable VPHP generators. Although this configuration is still com-

mon for smaller facilities and research isolators, many large production isolators now utilize integrated generators that work in concert with their air-handling systems to distribute and purge hydrogen peroxide, resulting in short, effective decontamination cycles. Once an isolator is cleaned and dried, hydrogen peroxide vapor is applied for biocontamination and to achieve sterility on all exposed hard surfaces.

Since the early 1990s, vapor-phase hydrogen peroxide use has expanded from barrier isolators to cleanroom suites. With the goal of achieving 4- to 6-log biocontrol, isolators have either purchased portable VPHP generators or contracted with specialized service providers that deploy multiple networked portable units, power cords, and fans to treat areas that can exceed 1000 cubic meters (m<sup>3</sup>). Gaseous biocontrol ventilation can be obtained more cost-effectively when VPHP is delivered via an integrated system that is seamlessly incorporated into the HVAC system of classified production areas.

For successfully integrated VPHP installation and optimal process performance, specialized user requirements may require a high level of cross-functional collaboration, with input from process engineers, automation programmers, quality, validation, and environmental health and safety personnel. Hundreds of integrated VPHP systems have been installed, primarily in Europe and Asia, providing clear evidence that this improved automated technology has gained industry acceptance.

### PLANNING CRITERIA

The decision to implement a VPHP process begins with selecting either a portable or integrated system. The choice often depends on existing or planned infrastructure, technical resources, intended frequency of use, and a favorable cost-of-ownership projection. A comparison is shown in Table A.

### Safety

Safety should always be a primary focus, even in the preliminary design phase. VPHP systems typically use concentrated liquid hydrogen peroxide in a range of 35%-65% strength and handling of these containers should comply with applicable rules. Planning and selecting VPHP from delivery, piping or an enclosure should also be a primary concern. Hydrogen peroxide breakdown into oxygen and water, accelerating the VPHP usually, is not an environmental concern. A risk assessment should be conducted to determine if workers or pedestrians near the exhaust outlet could be exposed to VPHP if necessary, catalytic converters can be installed to degrade VPHP before exhaust to outside the risk.

## Speaker

- > Working with VPHP applications at Steris Corp since 2004
- > Observations vs. formal testing
- > Difference between integrated and portable VPHP applications



Connecting


Pharmaceutical

Knowledge

ispe.org | 3

## Why Use VPHP?

- ✓ Efficacy (Broad spectrum sterilant)
  - Reach difficult to access surfaces
  - Passes through HEPA filters
  - Kills airborne and surface microbes
- ✓ Consistency & Distribution
  - Reach difficult to access surfaces
  - Passes through HEPA filters
  - Kills airborne and surface microbes
- ✓ Excellent Material Compatibility
  - Electronics
  - Metals and common polymers
- ✓ Speed
  - Minimal labor required
  - Easy to validate
- ✓ Green Technology
  - Low toxicity
  - No residues
  - EPA approved



Connecting

Pharmaceutical

Knowledge


ispe.org |

Agent	Uses
Chlorine Dioxide	Containment , industrial, water
Atomized Hydrogen Peroxide	Service decon
Vapor Hydrogen Peroxide	Service decon, cleanrooms, isolator, packaging
Ozone	Mostly food and water
Peracetic acid	Rooms, small enclosures
Ethylene Oxide	Mostly contracted service bulk product
Methyl Bromide	Agriculture , some interest from Homeland
UV light	Room decon, packaging
E-beam	Syringe tub decon

## 3 Key Factors for Robust and Rapid Cycles

Factor	Desired Condition
Percent Saturation	High percent saturation but below dew point
Concentration	Highest concentration without condensing
Distribution	Homogenous distribution

## Dose x Time = Kill



A D-value is the time needed to achieve a log reduction (90% kill) of the target organism




Connecting Pharmaceutical Knowledge

ispe.org | 7

## Efficacy More than Doubles with Percent Saturation (note: a single concentration used)

D-value for 1.6 mg/L versus Percent Saturation




Connecting Pharmaceutical Knowledge


ispe.org |

## Distribution - Chemical Indicators

### Less Distribution/Exposure



### More Distribution/Exposure



Connecting

Pharmaceutical

Knowledge

ispe.org | 9


## Most Common Applications



Portable  
Generator (piped)



Portable  
Generator (non-piped)



Integrated  
Generator (with manifold)






Enclosure



Airflow

### Scenarios

- Inside single or networked
- Outside – closed loop or single pass
- Integrated single pass piped to decontaminate zone 1 or zone 2



Connecting

Pharmaceutical

Knowledge

ispe.org |

## Which VPHP System?

### portable



- ✓ Spaces not yet defined
- ✓ Uses in different buildings
- ✓ Typically less than 6,000ft<sup>3</sup>
- ✓ Cycle time not a constraint
- ✓ Use of fans not an issue
- ✓ Less frequent use

### integrated



- ✓ Large and small spaces up to ~40,000 ft<sup>3</sup> / cycle + unit
- ✓ Same enclosures repeatedly
- ✓ Frequent use (chamber)
- ✓ Short cycle times
- ✓ Automated sequenced decontamination of multiple rooms



Connecting

Pharmaceutical

Knowledge

ispe.org | 11

## Common Applications

- **BSL2 / BSL3**
- **Grade A and B / ISO 5 and 7 Rooms**
- **Decontamination Chambers**
- **RABS**
- **Incubators**
- **Isolators**
- **Rooms containing LAFWBs or BSCs**




Connecting

Pharmaceutical

Knowledge

ispe.org | 12

## VPHP Modular System




Connecting

Pharmaceutical

Knowledge

ispe.org |

## Distribution Schematic




Connecting

Pharmaceutical

Knowledge


ispe.org |

## Control via Building Automation System Automatización





Connecting Pharmaceutical Knowledge

ispe.org |




Connecting Pharmaceutical Knowledge

ispe.org | 16



## Duct Interface




Connecting

Pharmaceutical

Knowledge

ispe.org | 19

## Manifolds



Manual valves



Actuated valves



Connecting

Pharmaceutical

Knowledge

ispe.org | 20

## Manifold



Connecting

Pharmaceutical

Knowledge

ispe.org |

## Biopharmaceutical Fermentation Suite

**Volume: 32,000ft<sup>3</sup> (900m<sup>3</sup>)**

**Ceiling height: 28ft (8,5m)**

**Single pass, No fans**

**6-log reduction**

**Cycle time 6 hours**

Cycle Phase	Time min.	Airflow	Injection g/min
Dehumidification	30	6 A.E./ hour	-
Condition	30	120 cfm	96
Decontamination	90	120 cfm	60
Aeration	210	40 A.E./ hour	-

A.E. = Air Exchange



Connecting

Pharmaceutical

Knowledge

ispe.org |

## Pass-Through Chambers / Transfer Hatches / Material Air Locks



Shown above



Enclosure v Volume ft <sup>3</sup>	Enclosure surface material	Injection rate Condition g/min.	Injection rate Decon g/min.	Decon time Min 6 log	Decon airflow ft <sup>3</sup> /min	PPM	Aeration airflow ft <sup>3</sup> /min	Total Cycle Time min.
460 (6x8x9.5'L)	Stainless	32	23	12	120	1000	765	45
175 (4x6x7'L)	Epoxy paint	12	9	8	40	950	1750	30



Connecting Pharmaceutical Knowledge

ispe.org |

## Decontamination by HVAC Zone



VHP M1000-T4

Proposed area of installation  
on AHU mechanical floor



Zone 1

Zone 2




Connecting Pharmaceutical Knowledge

ispe.org |

## RABS – (Restricted Access Barrier Systems)

Modular VHP systems can rapidly decontaminate both Active & Passive RABS and the rooms they are housed in




Source: Pharmaceutical International



Connecting Pharmaceutical Knowledge


ispe.org | 25



Connecting Pharmaceutical Knowledge

ispe.org | 26

## Integrated installations




Connecting

Pharmaceutical

Knowledge

ispe.org | 27




Connecting

Pharmaceutical

Knowledge

ispe.org | 28

## Cycle Development / Validation




Connecting

Pharmaceutical

Knowledge

ispe.org |

## Stakeholders - Define & Develop Your Expectations



Connecting

Pharmaceutical

Knowledge

ispe.org | 30

## Take Aways

### **VPHP –**

- **Can be deployed as a utility**
- **Automated decontaminated sequences can run from the BAS**
- **Once installed and validated – minimal effort to execute cycles**
- **Many different enclosure types can be decontaminated**  
isolators, RABS, cleanrooms, incubators, BSCs, Chambers  
etc.

**Once installed very low operating costs.**



Connecting

Pharmaceutical

Knowledge

ispe.org |

## Frequently Asked Questions

- **Does VPHP affect ductwork or room finishes?**
- **How is vapor contained in the target area?**
- **Cycle duration?**
- **Differences between an integrated and portable system?**
- **Enclosure qualification criteria?**
- **Most common mistakes?**





Connecting

Pharmaceutical

Knowledge

ispe.org | 32

## Questions ?



[John.klostermyer@steris.com](mailto:John.klostermyer@steris.com)

(908) 507-1131



Connecting

Pharmaceutical

Knowledge  
33

ispe.org |